



# Bio Plas Universal PCR Plates

## Bio Plas Universal PCR Plates

- Broad thermal cycler compatibility.
- High-contrast, black alphanumeric indexing for easy well identification.
- White wells enhance fluorescence reflection.
- Ultra-thin, uniform wells ensure optimal heat transfer for high reaction efficiency.
- Well ridges reduce sample evaporation when sealing with film, foil, or strip caps.
- Manufactured from high-clarity, ultra-pure, polypropylene to ensure purity and provide excellent optical clarity.
- Certified free of detectable DNase, RNase, DNA, PCR inhibitors, and tested pyrogen-free.

Cat# 9400



96 Wells. A12 Corner.

### Universal PCR Plate Ordering Information

Cat#	Description	Color	Pack Qty	Packs/Case
9400	0.2 mL Universal PCR Plate, Non-Skirted, Standard	Clear	10	10
9401	0.2 mL Universal PCR Plate, Non-Skirted, Standard	Blue	10	10
9402	0.2 mL Universal PCR Plate, Non-Skirted, Standard	Green	10	10
9403	0.2 mL Universal PCR Plate, Non-Skirted, Standard	Red	10	10
9404	0.2 mL Universal PCR Plate, Non-Skirted, Standard	Violet	10	10
9405	0.2 mL Universal PCR Plate, Non-Skirted, Standard	Yellow	10	10
9406	0.2 mL Universal PCR Plate, Non-Skirted, Standard	Asst	10	10
9407	0.2 mL Universal PCR Plate, Non-Skirted, Standard	White	10	10
9408	0.1 mL Universal PCR Plate, Non-Skirted, Low Profile	Clear	20	5
9409	0.1 mL Universal PCR Plate, Non-Skirted, Low Profile	White	20	5

Cat# 9408



96 Wells. H12 Corner.

Cat# 9509



### BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS

4380 Redwood Hwy, Suite C-16

San Rafael, CA 94903 USA

Tel: (415) 472-3777 Fax: (415) 472-3758

info@bioplas.com www.bioplas.com