



## Screw Caps For Microcentrifuge Tubes Material Safety Data Sheet

**Product Name:** Screw Caps For Microcentrifuge Tubes

**Bio Plas Catalog Numbers :** 4215 through 4224, 4215R through 4224R, and 4215RSB through 4224RSB.

**Company Name:** Bio Plas, Inc.

**Address:** 4380 Redwood Highway, Suite C-16, San Rafael, CA 94903

**Tel:** (415) 472-3777

**Recommended Use:** Capping screw cap microcentrifuge tubes.

**Product Composition:** Polypropylene. Free of RNase, DNase, DNA, Pyrogens and PCR inhibitors.

**Hazards Identification:** Not classified as dangerous for use/supply

**Hazardous Ingredients:** None

**Fire Hazard:** Product burns at 390 Degrees Celsius. The fumes are CO and CO2

**Fire Extinguishing Media:** Water, Foam, Carbon Dioxide, Dry Chemical

**Toxicology Information:** Material is non-toxic

**Storage Requirements:**

- This product may react with strong oxidising agents & should not be stored near such materials
- Storage temperature should be ambient. (preferably below 50 degrees C)
- Open flames prohibited.
- Protect from heat and direct sunlight.

**Note:** Bio Plas Screw Caps are only guaranteed to fit Bio Plas Screw Cap Tubes.