



Bio Plas Pestles

Bio Plas disposable polypropylene Pestles may be used to homogenize cells or re-suspend precipitates such as protein or DNA pellets in buffer solutions.

All the packaging options contain products free of detectable RNase, DNase, DNA and Pyrogens. These pestles and tubes are precision manufactured and quality controlled to firmly match together for maximum performance efficiency.

Pestles are conveniently packaged into the appropriate size tubes, ready for use and are sterilized.



1.5 mL Blue, 4030-PB

Ordering Information			
Catalog #	Description	Pkg. Qty	Pkgs./Case
4030-PB	1.5 mL Blue	10 per Bag	25 Bags
4040-PB	0.6 mL Blue	10 per Bag	25 Bags

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS
 4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA
 Tel: (415) 472-3777 Fax: (415) 472-3758
 info@bioplas.com www.bioplas.com