



BIO PLAS INC.
INNOVATIVE IDEAS...
FOR TODAY'S LABORATORY DISPOSABLES
MADE IN THE USA



Uni-Flex® Safety Caps

A safe, economical way to cap and recap blood collection tubes and disposable glass culture tubes, regardless of manufacturer. **Uni-Flex® Safety Caps** are available for a wide range of biomedical and industrial applications.

The **Uni-Flex® Safety Cap** liquid tight seal ensures sample integrity and prevents spillage of serum samples. **Uni-Flex® Safety Caps** protect your laboratory staff from the aerosols of highly infectious organisms such as TB and the HTLV-III virus.

Uni-Flex® Safety Caps Benefits:

- Will not aerosol when opened
- 10, 12, 13 and 16mm sizes fit all varieties of blood collecting and disposable culture tubes
- Liquid tight seal prevents sample evaporation. The **Uni-Flex® Safety Cap** is designed to hold under centrifugation, agitation, refrigeration and freezing.
- Easy identification of samples: Uni-Flex® Safety Caps are available in seven colors that conveniently correspond to blood collection cap colors
- The **Uni-Flex Safety Cap** locks securely onto the outside of the tube lip to prevent evaporation.
- Easy, one hand removal and application. To close, press down firmly with thumb to lock in place. To open, simply pry cap up with thumb.
- Easy to label: Write directly on the cap with indelible ink
- Clear, natural colored caps are trace-metal free



Three Internal Seal Rings



BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS

4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA

Tel: (415) 472-3777 Fax: (415) 472-3758

info@bioplas.com www.bioplas.com

Uni-Flex® is a registered trademark of Bio Plas, Inc., San Rafael, CA



BIO PLAS INC.
 INNOVATIVE IDEAS...
 FOR TODAY'S LABORATORY DISPOSABLES
 MADE IN THE USA



Ordering Information

Catalog #	Description	Color	Pkg. Qty.	Pkgs./ Case
6500	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	White	1,000	10
6505	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Red	1,000	10
6510	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Lavender	1,000	10
6515	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Green	1,000	10
6520	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Blue	1,000	10
6525	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Grey	1,000	10
6530	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Yellow	1,000	10
6550	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	White	1,000	10
6555	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Red	1,000	10
6560	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Lavender	1,000	10
6565	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Green	1,000	10
6570	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Blue	1,000	10
6575	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Grey	1,000	10
6580	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Yellow	1,000	10
6600	Uni-Flex®Safety Caps for 13mm Culture Tubes	White	1,000	10
6605	Uni-Flex®Safety Caps for 13mm Culture Tubes	Red	1,000	10
6610	Uni-Flex®Safety Caps for 13mm Culture Tubes	Lavender	1,000	10
6615	Uni-Flex®Safety Caps for 13mm Culture Tubes	Green	1,000	10
6620	Uni-Flex®Safety Caps for 13mm Culture Tubes	Blue	1,000	10
6625	Uni-Flex®Safety Caps for 13mm Culture Tubes	Grey	1,000	10
6630	Uni-Flex®Safety Caps for 13mm Culture Tubes	Yellow	1,000	10
6700	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	White	1,000	10
6705	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Red	1,000	10
6710	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Lavender	1,000	10
6715	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Green	1,000	10
6720	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Blue	1,000	10
6725	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Grey	1,000	10
6730	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Yellow	1,000	10



Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes, Cat# 6525



Uni-Flex®Safety Caps for 13mm Culture Tubes, Cat# 6620



Uni-Flex®Safety Caps for 13mm Culture Tubes, Cat# 6605

“The Premium Tube Cap With The Triple-Lipped Internal Seal!”

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS

4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA

Tel: (415) 472-3777 Fax: (415) 472-3758

info@bioplas.com www.bioplas.com

Uni-Flex® is a registered trademark of Bio Plas, Inc., San Rafael, CA